Asphaline temp

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Asphaline temp

Composition:

Polymerically adsorbed iodine, Zinc oxide, Barium sulphate

The paste contains 3% by weight polymerically bonded iodine.

Properties:

Medical product for the temporary filling of gangrenous root canals as an auto-sterile barrier for protection against re-infection from the oral cavity and to improve the healing process before permanent sealing. Non-toxic and non-irritant. Odourless with neutral taste. No interactions or resorption by the tissue are to be expected. The paste has very good shear adhesion. No pore formation. Easy working of the ready-to-use paste and retreatment using the familiar techniques of root canal treatment.

How to use:

Prepare cavity and root canal, rinse and dry. Rinse out any residual hydrogen peroxide and hypochlorite solution, since it may slightly impair the effect. Select lentulo and coat with a little paste. Work the paste around slowly into the root canal. Repeat procedure if required. Do not overfill. You can check by X-ray that the filling is complete.

Asphaline temp does not set completely, and should be removed after a period of 5 days. With repeated use, an accumulated period of 14 days in the root canal should not be exceeded.

Restrictions on use:

Asphaline temp should not be used in patients that have an allergic reaction to one of the product constituents or have had a thyroid disease or take lithium-based psychopharmaca. It must not be used during pregnancy and lactation.

Instructions and guidance:

Do not warm paste up. The product must only be used by dentists. Stains on the skin and clothing generally disappear after 24 hours, although they can be removed using sodium thio-sulphate solution.

If it comes into contact with blood no reactions or limitation on effectiveness are to be expected. Iodine decomposes into iodide and iodate when dissolved in aqueous solution. Irritation of the periapical tissue are therefore not to be expected. The colour of the material may fade away by the discharge of iodine.

The total quantity of iodine used for each root canal (approx. 5 µl in volume) is 150 µg. The proportion that is discharged whilst the paste is in the root canal falls well below the recommended daily dose of 100 – 200 µg.

Keep Asphaline temp in a cool and dry place. Do not expose to direct sunlight. Close properly after use. After opening the paste may show a whitish or yellowish colour at the tip of the tube. The discoloured material should be removed. Once started, use tube up within 6 months. Do not use after the expiry date stated on tube and packaging. Keep out of the reach of children.

The paste can be mixed together with aqueous materials.

Always quote the lot number whenever the material has to be identified.

REF 550

Tube with 7 g ready-to-use paste

The product was developed for use in dentistry and must be used only in accordance with the instructions for use. The manufacturer accepts no liability for loss or damage caused by improper or inappropriate use. The possessor himself has to verify whether the product is suitable for the use he intends, should there be no corresponding information in the instructions for use.

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Alfred Becht GmbH